

# SomnoGuard®

AP2 & SPX

Mandibular Advancement Devices

for OSA  
& Snoring

## Effective, Efficient, Affordable Oral Appliance Solution

- Designed to be an integral part of **YOUR** practice
- Custom Fit in a single 15-20 minute office visit
- Full Mandibular and Maxillary Support
- Precisely Titratable to optimize efficacy
- Satisfies Medical Necessary criteria of many insurers
- Vastly more affordable - even when not covered

**#1** Therapy solution for Snoring

**#2** Therapy solution for Obstructive Sleep Apnea  
- 2nd only to CPAP



Precision Manufactured in Germany by

**TOMED**  
medical products

**1ST**  
**LINE** Medical Inc.

P : 866-720-8080

E : Service@1stLineMedical.com

W : www.1stLineMedical.com

O : 850 Route 3, Ste 201

Holderness, NH 03245

US Importer/Distributor

**74%**

COMPLIANCE

**65%**

AHI  
REDUCTION

# PROVIDE YOUR OWN SOLUTION

Improve patient care with:

Greater Access

Expedited Care

Lower Healthcare Costs

Maintain Patient Relationship

Build Sleep Practice

Continuum of Care

## Custom Fit Prefabricated Appliances



### SomnoGuard AP2

PN 03050200

Low profile design

2-part, 2-piece, one size fits most

Freedom of lateral movement optimizes  
comfort/compliance

Micro-Titratable: 0-10mm



### SomnoGuard SPX

PN 03090000

Ultra-low profile

2-part, 2-piece, variable width accommodating  
enhanced retention copolymer

titration mechanism permits mouth breathing

Micro-Titrating turnbuckles or 1mm increment straps


SomnoGuard<sup>®</sup>  
MANDIBULAR  
ADVANCEMENT DEVICES  
AP2 & SPX

# E0485

HCPCS CODE

# G47.33

ICD-10 CODE



Scan for  
video presentation

## Dispense or Prescribe - it's up to you

1. **Purchase** (\$199-\$259) and dispense (MSRP \$700)
2. **Prescribe** (\$259-\$279 +\$8.95 S/H) Patient Purchases and schedules billable office visit for custom fitting (ie. CPT 99213)
3. Or **both**, depending upon the circumstances